

ACT ONE

# The amazing power of scripting DS9, and some totally rad things you can do with it.

Establishing the station.

## Cooler than it sounds, really!

INT. JAKE & NOG'S QUARTERS

Nog is sitting in a chair with his arms crossed,  
glaring at Jake and determined not to be swayed. Jake  
is just as determined to get the diminutive Ferengi  
Cadet to bend.

NOG

No.

JAKE

Come on, Nog...  
By


NOG



No.

Kyle Kaplan  
Why not?




# What is DS9?

<http://ds9.si.edu/site/Home.html>







-20 -4 12 30 48 66 84 102 120




0.8 2 0.28 1.4 3.2 4.7 6.8 8.7



4.4 6.8 12 18 24 30 36 40



# Python



# XPA

- <http://hea-www.harvard.edu/RD/xpa/>
- Allows you to control DS9 via the terminal
- Open DS9 from command line
  - *ds9 &*
- Send command to DS9:
  - *xpaset -p ds9 blahblahblah*
- Grab information from DS9:
  - *xpaget ds9 blahblahblah*
- More info on XPA commands: <http://ds9.si.edu/ref/xpa.html>

# Python Library Wrapper for XPA & DS9

```
#Python library for accessing DS9 and XPA.  
#Written by Kyle Kaplan March 2014.  
from subprocess import call #Allow python to access command line  
from subprocess import check_output #Allow python to access command line and return result to a variable  
import time #To put in delays  
  
#Open DS9  
def open():  
    call('ds9 &', shell=True)  
    time.sleep(2) #Give computer a few seconds to respond after opening DS9  
  
def close():  
    call('xpaset -p ds9 exit')  
  
#Get xpaget statements from ds9  
def get(command):  
    result = check_output('xpaget ds9 ' + command, shell=True) #Get information from ds9 using XPA get  
    return str(result).strip() #return information grabbed from ds9  
  
#Send xpaset commands to ds9  
def set(command):  
    call('xpaset -p ds9 ' + command, shell=True)  
    #time.sleep(0.5) #Set delay after each command to give computer time to respond before the next command  
  
#Allow user to set delays in DS9 scripts  
def wait(delay):  
    time.sleep(delay)  
  
#Special command for drawing regions onto ds9  
def draw(command):  
    print 'echo "' + command + '" | xpaset ds9 regions'  
    call('echo "' + command + '" | xpaset ds9 regions', shell=True)  
  
def rot(angle):  
    set('rotate '+str(angle))  
  
def rotto(angle):  
    set('rotate to '+str(angle))  
  
def north():  
    set('rotate to 0')
```

# Demo Time!

- Demos
  - Rotate
  - Zoom
  - Compare catalogs
  - Grabbing and displaying crosshair coordinates
  - Make a DSS finder chart
- Real world scripts I have made
  - Find brightest stars in field
  - IGRINS finder charts